

Better for the Baby

Better for the Nurse

Better For The NICU/Hospital

Better X 3

























www.small-beginnings.com 1-800-676-0462 Fax orders to 1-760-948-1916



This page is dedicated to Nurses, who have supported us and helped us bring our products to the neonatal intensive care units around the world. Your desire to do your best for the infants in your care has been the reason Small Beginnings Inc. is considered the leader in producing great patented neonatal devices. When you called us with needs or problems about products that didn't work, we have rolled up our sleeves and created solutions.

Our Mission Statement is as true today as it was when we started our company.

We will make products that are better for the Baby! We will make products that are better for the Nurse! We will make products that are cost-effective and better for the hospitals! Better x 3.

It is this philosophy that has individual hospitals and entire hospital systems converting to our wonderful product line. We are proud to bring you the very best year after year and with your support we will continue to do so.

My Sincere Thanks

M. Sharon Rogone RNC-E Chief Executive Officer Small Beginnings Inc. A Certified Women



**Owned Business** 

	rable of	Contents		
category 5 Item	Pg#	category 5 Item	P9#	
Fantastic Four		*continued		
Bili-Bonnet	4	Dríp-drop Solutions		
Pacee Pacifier	24	Bebe-Chix	34	
Bebeonkers	26	Birth N Bath Wraps/Bags	35	
Diapers	36, 38	Cuddle Buns Diapers	36	
Suctioning Solutions	76	Beary Bili Buns Diapers	38	
Bebeonkers 0, 1, 11, 111	26	Other solutions		
Bebe-Ball	28	Climate Cover	6	
Bebe-Keeper Kover	29	Quillo-Iso Cover	2	
Positioning Solutions		Pacee Pacifier	24	
Heavenly Mattress	10	Parent Guides	30	
Wonder Star	12	Preemie Gowns	31	
Cuddle-Up	14	Wing-slings	32	
Nests	16	Privacy Poncho	33	
Bed Buddy	18	New Solutions		
Posture Pillows	20	Neck rolls	40	
Píllow Sampler	22	Birth N Bath Wraps/Bags	35	
Positioning Sampler	23			

#### Bili-Bonnet





Stretches open to allow Bili-light to penetrate the entire head.

3 Sizes Fit all head shapes.

- No Confusion
- Easy to apply
- No straps to frustrate
- No adhesives to damage skin



Compression Technology

Makes ours the only mask with
minimal light penetration

- Rests on ocular face bones and removes pressure from eye
- Mask is Reversible made of PE Foam & Brushed Nylon
- Conforms to baby's face with plenty of room for adjustment on each side

Bili-Bonnet Phototherapy Masks

O1A Preemie size for infants approx. 27-30 weeks

O1B Newborn size for infants above 31 weeks approx.

Ultra preemie size for infants below 27 weeks approx.

Tabs are rounded & very low profile hook & loop that allow attachment at any angle to maximize fit and minimize slippage.

## Bili-Bonnet

# Latex Free-No adhesives

#### How to put the patented Bili-Bonnet phototherapy mask on correctly



Remove from package, then stretch gauze and pull over fist to create a bonnet.



Stretch bonnet over baby's head to go from the nape of the neck to the forehead.



Attach velcro tab to one side.

Note, the mask has a 1 inch
area for better adjustment &
velcro tabs are rounded.



Repeat with tab on other side. The mask should be snug & secure but not tight.



Ensure the netting is opened up all the way over the head to allow light to penetrate.



Now once again stretch netting down to the nape of the neck and up to the forehead.



Most Important Step
Once all is in place and
secure, unhook and reattach
the velcro tabs to pull
up/down/midline as head
shape dictates.



The end result is a mask that stays on better giving the nurse more time for other things.

Protects the eyes better & allows more light to the head.

Take time to save time!

## Climate Cover III




"The perfect micro-environment"

- Reduces insensible water loss and cold stress
- Reduces energy cost as it maintains body temperature
- Becomes an OXY-BODYHOOD when infant is off the ventilator
- \*\*Warning must use with O2\*\*

Order #	Product Description	B× Qty
	Climate Cover III	
	The imoregulation	
02E	Oxy-Body Hood	6
	19×25×10	

- · Can stay in place for most care with easy access Arm Ports
- Nurses can do good developmental care and positioning
- · Can be X-rayed through
- Maintains a neutral thermal environment even if/when heat source must be removed for a brief time






For removal in case of emergencies or other procedures all lines and tubing go under the bottom arches on all 4 sides

## Climate Cover III





Can make every radiant warmer a

"Poor-man's Giraffe"



## Iso-Sphere "Quillo"

Universal Isolette Cover

Infants in the NICU are bombarded with constant noise, harsh lights and painful procedures. These tiny babies learn quickly that when they are covered they do not get poked, prodded and stuck. Reducing stress promotes faster healing, ideal growth and proper development which in turn reduces hospital stays with better outcomes.

- 8 ounce quilting to reduce noise & light
- Machine washable and reusable
- Cotton/poly blend & flame resistant
- Gray non-stimulating inside reduces stress and promotes better rest/sleep These covers serve a great purpose!!





Order # Product Description

By Otty

Iso-Sphere Cover

03DD

Universal Isolette Cover Fold into a pocket "Quillo"

6

These flame resistant covers not only fit over any isolette, bassinette or even cribs, they also easily cover the "Giraffe" as well. When not in use it folds into itself to make a pillow or "Quillo". It can also be labeled with your hospital and/or the baby's name and birth date for a perfect keepsake from the NICU.

Can also be carried in your hospital's gift shop.



## Heavenly Foam





- Visco-Elastic, first developed by NASA
- Latex free memory foam
- · Perfect topper for the neonatal bed
- Reduces surface tension, pressure points & stress
- In 2 sizes and is 5/8 of an inch thick



2 Straps to Use as Needed

> Heavenly Foam is included and used with Cuddle Up. Note use with other products as well as Bed Buddy to provide needed soft rebound and boundaries.

Recommend Use with Bed Buddy





Vinyl Covered Heavenly Foam with Knit Cover

Heavenly Foam Mattress			
Order #	Product Description	Cs Qty	
	Heavenly Foam Mattress only	"	
04VE	Preemie/Newborn Heavenly Foom Mattress - 20"×13.5"×5/6"	24	
04CVE	Ultra Preemie Heavenly Foom Mattress 10"×13.5"×5/8"	10	
	Leatherette Vinyl Covered Heavenly Foam		
04VE-C	Preemie/Newborn Vinyl Covered mattress	В	
04CVE-C	Ultra Preemie Vinyl Covered mattress	6	
04A VE-C	Vinyl Covered mattress for Preemie/Newborn with Fleece Covers	8	
04AKVE-C	Vinyl Covered mattress for Preemie/Newborn with Knit Covers	В	
04CCVE-C	Vinyl Covered mattress for Ultra Preemie with Fleece Covers	6	
04CCKVE-C	Vinyl Covered mattress for Ultra Preemie with Knit Covers	6	

"We only make products that are better for the baby, better for the nurse and better for the hospital." – M. Sharon Rogone RNC-E CEO

## Heavenly Foam Covers

- · Various coverings available
- Heece, knit & disposable
- Waterproof disposable
- Vinyl covered Heavenly Foam also available





Waterproof Disposable Cover



Order #	Product Description	Cs Qty
	Covers for Vinyl Covered Heavenly Foam Mattress	
04A	Preemie/Newborn Fleece Covers	10
04AK	Preemie/Newborn Knit Covers	10
04C	Ultra Preemie Pleece Covers	10
0.4CK	Ultra Preemie Knit Covers	10
04VED-C	Preemie/Newborn Disposable Covers	20/Pkg
04CVED	Ultra Preemie Disposable Covers	20/Pkg
	Disposable Covers for Heavenly Foam Mattress	
04VED	Disposable covers for Preemie/Newborn Mattress	20/Pkg
04CVED	Disposable covers for Ultra Preemie Mattress	20/Pkg
	Baby Blankee	
04b	Boby Blorket 25"x31"	12

## Wonder Star

#### Wonder Star

- Multi-purpose positioning device
- Soft gel rebound beads
- Washable
- Move & manipulate in many ways
- Keep baby comfortable.
- 2 sizes
- Fleece or Knit
- Great addition to any NICU.



Order #	Product Description	Cs Qty
Wonde	er Ster Multi-Purpose Positioning Pillo	)W
05A	Wonder Star in Fleece 15"arm span and 8" wide body	10
05AK	Wonder Star in Knit 15"arm span and 8" wide body	10
050	Lil' Wonder Star in Fleece 13"arm span and 6.5" wide body	10
05 <i>C</i> K	Lil' Wonder Star in Knit 13"arm span and 6.5" wide body	10



## Wonder Star





Knit Star

Very flexible positioning

Available in Fleece or Knit this
multi-functional positioning aid can meet
many needs. Providing the best positioning
for your baby, whether maintaining
head or body position or nestling baby in a
womb-like position.

The Wonder Star can do it all.






## Cuddle Up

The Wonderful Cuddle Up
Comes with a Heavenly Foam
mattress to put in its
pillow case sleeve.
Recommend use with the
Bed Buddy to provide soft,
push back boundaries

Nurses requested a better swaddling device, so Sharon Rogone created the Cuddle Up Swaddler.

Order #	Product Description	Cs Qty
	Cuddle Up Sets	
DBA	Preemie/Newborn Fleece Guddle Up with Foam Mattress	6
DBAK	Preemie/Newborn Knit Cuddle Up with Foam Mattress	6
Dec	Utira Preemie Guddle Up with Foam Mattress and 30°Bed Buddy	LO
DECK	Utira Preemie Gudd le Up with Faam Mattress and 30° Bed Buddy	LO
	Cuddle Up Enugglers	
065 A	Preemie Cuddle Up Snuggler	6
065 B	Newbarn Guddle Up Snuggler	6
065 C	Uttra Preemie Cuddle Up Snuggler	6
065 CW	Ultra Preemie Cudd le Up Snuggler Warming Winap	6

Introducing the new Snuggler wrap
designed to maintain infants
temperature and secure tubing and
lines during any procedure; such as,
transports, weighing, exams, and
kangaroo care, etc. Is also available
in ultra preemie size only, with
2 pockets on the back side for
placement of heel warmers.





## Cuddle Up Covers

Fleece Cuddle Up Cover



Note 2 straps

With inner sleeve for the Heavenly Foam Mattress.



Disposable Cuddle Up Cover

Order #	Product Description	Ćs Qty
	Cuddle Up Covers only	
06AA	Preemie/Newborn Fleece	10
06AAK	Preemie/Newborn Knit	10
06CC	Ultra Preemie Fleece	10
06CCK	Ultra Presmis Knit	10
06AAD	Preemie/Newborn Disposable Cover	20/Pkg
0.6CCD	Ultra Preemie Disposable Cover	20/Pkg

Cuddle Up Covers

- Fleece, Knit or Disposable Covers
- 2 straps to Cuddle, snuggle and secure
- · Sleeve for Heavenly Foam Mattress

How do you keep the Heavenly Foam mattress pad from being soiled or wet?
The answer is our waterproof disposable covers lined with plastic to fix the problem.
Now your baby can lay on NASA's best for comfort and circulation
and know it is clean and dry.

Waterproof Disposable Cover





#### Nest Sets

#### Created for the sensitive, hard to manage infant

- Decrease decels and desats
- Reduce stress and negative reactions & set backs
- Gel beads push back for spatial & brain development
- Bags support baby's back, front or side
- Provides for proper positioning
- 3 sizes: Preemie, Newborn and Ultra Preemie to meet the needs of all age & size babies



Super soft gel beads

Patented by our founder, the Preemic Nest replaces rolled blankets & time lost by the nurse



The ultimate swaddler
& positioning device

We added a 4th sleeve so you can rearrange the bags & only turn baby

Now with Velcro tabs to keep bags from sliding out



Order #	Product Description	Cs Qty
	Nest Sets Infant Positioning Device	
09C	Preemie Nest Set with Fleece Pad	10
09CK	Preemie Nest Set with Knit Pad	10
09CC	Newborn Nest Set with Fleece Pad	6
09CCK	Newborn Nest Set with Knit Pad	6
09CCC	Ultra Preemie Nest Set with Fleece Pad	10
09CCCK	Ultra Preemie Nest Set with Knit Pad	10

#### How to prepare the Nest

- Center baby in supine position
- Next grab bead bag close to you
- Roll the bag toward yourself to turn baby from supine to side position
- Push other bead bags in against baby
- Surrounds baby with womb-like feeling, soft, cushy, rebound
   Turning & Positioning in seconds!!









The Nest Pads are designed with ribbed knit sleeves sewn into the back of a super soft fleece or knit blanket. Each sleeve then has Velcro at each end to keep the bags from falling out. Pads come with two long and two short sleeves for ease in turning and positioning. Disposable pads are also available



Veloro fastening to hold bead bags in place

Order #	Product Description	Cs Qty
	Nest Pads only	
09A	Preemie Fleece Pods	12
09AK	Preemie Knit Pods	12
09AA	Newborn Fleece Pads	12
09AAK	Newborn Knit Pads	12
09AAA	Ultra Preemie Fleece Pads	12
09AAAK	Ultra Preemie Knit Pads	12
7.1	Disposable Pads Only	Pkg Qty
09AAD	Preemie/Newborn Disposable Pads	20
09AAAD	Ult ra Preemie Disposable Pads	20

## Bed Buddy

The Bed Buddy is the easiest positioning pillow to use. It quickly surrounds and interfaces with infant, providing a soft rebounding surface. It has zero replacement cost. Just wipe vinyl clean with any disinfectant, wash and replace cover and it's






Wipe clean vinyl bead bag

	Bed Buddy Infant Positioning Devices		
Order #	Product Description	Cs Qty	
	Bed Buddy with Covers		
144	Presmis 44" Bag with Fleece cover	6	
14AK	Preemie 44" Bog with Knit cover	6	
148	Newborn 60" Bag with Fleece cover	6	
148k	Newborn 60"Bag with Knit cover	6	
14C	Ultra Preemie 36" Bag with Fleece cover	10	
14Ck	Ultra Preemie 36" Bag with Knit cover	10	
24E	XS Ultra Preemie 30"Bag with Pleece cover	10	
14EK	XS Ultra Preemie 30" Bag with Knit cover	10	
50	Bed Buddy Bag only		
14.4-8	Preemie 44"Bag	6	
148-8	Newborn 60"Bag	6	
14C-B	Ultra Presmis 36" Bag	10	
14E-8	XS Ultra Preemie 30" Bag	10	

## Bed Buddy Covers




Velcro tabs to secure around infant

Fleece Cover

Order #	Product Description	Cs Qty
	Bed Buddy covers only	
1444	Presmis Flasce Covers	10
14AAK	Preemie Knit Covers	10
1488	Newborn Fleece Covers	10
1488K	Newborn Knit Covers	10
1400	Ultra Preemie Fleece Covers	10
14CCK	Ultra Preemie Knit Covers	10
14EE	XS Ultra Preemie Fleece Covers	10
14EEK	XS Ultra Preemie Knit Covers	10
	Disposable Covers only	Pkg Qty
144 AD	Preemie Disposable Covers	20
1488b	Newborn Disposable Covers	20
14CCD	Ultra Preemie Disposable Covers	50
14E ED	XS Ultra Preemie Disposable Covers	20



The Bed Buddy can be covered with our washable Fleece or Knit and is also available in our print disposable.



Disposable Cover

## Posture Pillows

Gal Baseds or Visco Bastic (VE)

Order #	Product Description	Cs Qty
	Form Posture Pillows with Covers	**
10A-VE	Preemie Abdominal Pillow with Fitted Fleece Cover	10
10AK-VE	Preemie Abdominal Pillow with Fitted Knit Cover	10
10AA-VE	Ultra Presmie Abdominal Pillow with Pitted Fleece Cover	10
10AAK-VE	Ultra Preemie Abdominal Pillow with Fitted Knit Cover	10
108-VE	Presmie Head Pillow with Fleece Cover	10
108K-VE	Preemie Head Pillowwith Knit Cover	10
1088-VE	Ultra Preemie Head Pillow with Fleece Cover	10
1088K-VE	Ultra Preemie Head Pillow with Knit Cover	10
108N8-VE	Newborn Head Pillow with Pleece Cover	10
108NBK-VE	Newborn Head Pillow with Knit Cover	10



#### Head Pillow-Gel Beads or VE Foam

- Allow better shaping of head
- Beads can simply be moved to one side to make a neck roll
- Thick VE memory foam you can x-ray through

Gel Bead Pillows



#### Abdominal Pillow-Gel Beads or VE Foam

- Softer TPE beads that feel like gel but never evaporate
- · High profile VE foam for greater height and can be x-rayed thru
- · Both developmentally designed with fitted covers
- Better prone positioning, better lung expansion & oxygenation

	Bood Bag Pasturo Pillows with Covers	
LDA	Preemie Abdamma   Pillaw with Fitted Fleece Caver	LD
LOAK	Preeme Abdammal Pillow with Frited Knit Caver	LD
LDAA	UhraPreemie AbdamnalPillowwith Fitted Fleece Cover	10
MAAK	Utra Preeme Abdomnal Pillow with Frited Knit Cover	LD
LDB	Preemie Head Pillow with Fleece Cover	LD
LOOK	Preemie Head Pillow with Knit Caver	LD
1088	Utira Preemie Head Pillow with Fleece Cover	LD
LOBBK	Ultra Preemie Head Pillow with Knit Cover	10
100	Preemie Wedge Pillow with Fleece Cover	10
LOCK	Preemie Wedge Pillow with Knit Cover	LD
1000	Ultra Preemie Wedge Pillow with Fleece Caver	LD
LDCCK	Ultra Preemie Wedge Pillaw with Knit Caver	LD



#### Wedge Pillows

- Multi-Functional
- Holds vent tubing & IV lines
- Quickly changes baby's position
- Wedge Pillow is a work horse
   You're sure to find your own uses!

## Posture Pillow Covers

#### Fleece \*\*

- Very soft
- Light-weight
- 100% polyester
- Machine washable
- Cheery colorful
- Multi-animal print

#### Knit \*\*

- Very soft
- 50/50 cotton/poly
- Machine washable
- Cheery colorful
- Multi-animal print
- White background

#### Disposable \*\*

- Soft
- Sturdy
- Non-woven
- Waterresistant
- Cheery colorful
- Multi-animal print

Order #	Product Description	Pkg Qty
	Posture Pillow Covers Only	
10E	Preemie Head & Wedge Pillow Fleece Covers	10
10EK	Preemie Head & Wedge Pillow Knit Covers	10
10EA	Presmie Abdominal Pillow Fitted Fleece Covers	10
10EAK	Preemie Abdominal Pillow Fitted Knit Covers	10
10EE	Ultra Preemie Head & Wedge Pillow Fleece Covers	10
10EEK	Ultra Preemie Head & Wedge Pillow Knit Covers	10
10EEA	Ultra Preemie Abdominal Pillow Fitted Fleece Covers	10
10EEAK	Ultra Presmis Abdominal Pillow Fitted Knit Covers	10



<ul> <li>Ab cover is fitted to new updated</li> </ul>	1 shape
---	---------

- Wedge coverdesigned to fit wedge pillow
- Head cover designed to fit head pillow



Order #	Product Description	Pkg Qty
	Disposable Covers	
10ED	Preemie Head & Wedge Pillow Disposable Covers	50
10EAD	Preemie Abdominal Pillow Disposable Covers	20
10EED	Ultra Preemie Head & Wedge Pillow Disposable Cove	20
10EEAD	Ultra Preemie Abdominal Pillow Disposable Covers	50



Packaged Posture Pillow Disposable Covers

## Posture Pillow Sampler Kits

\*\*\* Reference Article
Premature Infants: A Later Look
NYU Child Study Center

"Positioning the preemie can also influence future development. For example, preemies may be routinely left to lie on their backs, a position that differs from the natural curved fetal position they would maintain in utero. This can result in a shortening of the scapula muscles between the shoulder blades, causing a restricted range of motion and later motor problems. Thus, the older infant may be unable to appropriately reach for objects in front of him..."



Posture Pillow Sets Includes 1 of each pillow:

- Abdominal/prone positioning
- Head/neck or head
- Wedge/all around
- . With Knit or Fleece covers



Knit Cover



Disposable Cover





1	7
1	
1	Chide

-		4	
2	Sizes	Avai	loo le

Company of the Company		2000 M. M. M.
	Posture Pillow Sets with Gel Beads	
10b	Set: 1 Preemie Ab, Head, & Wedge Pillow & Fleece Covers	10 Sets
10DK	Set: 1 Preemie Ab, Head, & Wedge Pillow & Knit Covers	10 Sets
10bb	Set: 1 Ultra Preemie Ab Head & Wedge Pillow & Fleece Covers	10 Sets
10DDK	Set: 1 Ultra Preemie Ab, Head & Wedge Pillow & Knit Covers	10 Sets





All the positioning aids
any unit could ever want or need.
Every situation met and enough
for more than one baby.
All bags wipe clean, are reusable
and come with washable or
disposable covers.
Guaranteed for 2 years.
Covers come in knit, fleece or
disposable material.

All the amazing benefits of our 7 positioning products.
Assists PT, OT & all, with positioning, alignment, comforting and swaddling.
Remember when we position, comfort and reduce stress we are enabling normal growth and development.
When we mimic the womb erwironment we help the cognitive, spatial and physical growth of our tiny charges.
Never forget our oath to "Do No Harm."

Order #	Product Description	Cs Qty
	Positioning Sampler Kits	
07A	Preemie Positioning Sampler Kit includes 1 of each device with fleece covers: 10A Abdominal Pillow, 10B Head Pillow, 10C Wedge Pillow, 09C Nest Set, 06A Cuddle-Up, 05A Wonder Star, 14A Bed Buddy.	1 each/ 1 set
07AK	Preemie Positioning Sampler Kit includes 1 of each device with knit covers: 10A Abdominal Pillow, 10B Head Pillow, 10C Wedge Pillow, 09C Nest Set, 06A Cuddle-Up, 05A Wonder Star, 14A Bed Buddy.	1 eoch/ 1≤et
076	Ultra Preemie Positioning Sampler Kit with 1 of each device with fleece covers: 10AA Abdominal Pillow, 10BB Head Pillow, 10CC Wedge Pillow, 09CCC Nest Set, 06C Cuddle-Up, 05C Wonder Star, 14E Bed Buddy.	1 each/ 1 set
07CK	Ultra Preemie Positioning Sampler Kit with 1 of each device with knit covers: 10AA Abdominal Pillow, 10BB Head Pillow, 10CC Wedge Pillow, 09CCC Nest Set, 06C Cuddle-Up, 05C Wonder Star, 14E Bed Buddy.	1 each/ 1 set

## Pacee Pacifier





- 100% Silicone all one piece as recommended by The Journal of Pediatrics.
- 6 Large ventilation holes makes it the safest pacifier at discharge and for home.
- 2 sizes in easy to see colors
- · Teal for newborns and purple for preemies.






Order #	Product Description	B× Qty
	Pacee Pacifiers	
15A	Lil' Pacee preemie pacifier 100% silicone	100
15B	Nu'Paces newborn pacifier 100% silicons	100

## Pacee Pacifier









Use shallow dip for CPAP

Use deep V for ET tube

- Never needs cutting
- · Use deep V for ET Tube
- Use shallow dip for Nasal CPAP
- Collapses easily to the base
- Mimics nutritive suck



- · Softer nipple for premature infants
- · Restful non-nutritive suck
- · Won't easily fatigue infant
- Recommended by Lactation Consultants



#### Bebeonkers

Bebeonkers II

Bebeonkers II (16-2) Our original patented design is an Oral/Nasal suction device that does not grab tissue and effectively removes the most tenacious secretions, quickly and without compromising the infant.

Packaged six per bag.





Flexible & Soft

Bebeonkers I (16-1) was designed for the Ultra Preemie infant to facilitate quick removal of secretions both nasally and orally but still retains a three-hole configuration for safe oral care without hurting baby.

		man date
Order #	Product Description	B× Qty
	Bebeonkers Onal/Nasal Suction Device Only	
16 1	Bebeankers [ Uhra Preemie Bu bed Tip 6 pcs/z p bag,50 bgs/cs	300
16 2	Bebeankers [[ Standard Tp 6 pcs/xp bag, 50 bgs/cs	300
16-1KK	Bebeankers C with kover 3 Bebeankers +3 cover/bag, 50 bgs/cs	300
14-288	Bahasakasa (Custh kaya a 3 Bahasaka sa + 3 cayas/haa AB har/ca	300

Bebeonkers I

#### Bebeonkers





Devices altered to facilitate the removal of oral secretions or designed ineffectively for tenacious secretions. Thus drive we remove the loopists used to every as Nares sale whatever they can to action any electricity.

Major liability

All above are single hale drive all good base when sectioning and that's ball.

# remember a single hale device all good base when sectioning and that's ball.

Patented
Zor more holes
Bebeonkers has 3
Does not grab tissue
No trauma
No harm

Because there has never been an appropriate device for suctioning intubated or newborn infants, doctors, nurses and therapists have been using and creating inferior one-hole devices. They endanger the infants and make the hospital and medical professional liable if an infant is injured.

The patented

"BEBEONKERS" is that right device!!

It is the right tool for the right job.

Two or more holes and a tapered inner wall does the job quickly and without tearing tissue.

"Sometimes you just can't stand by when things are being done that can hurt a baby. We saw a problem that was international and made a product to solve it!"

—Philip N. Rogone PA. RCP.

#### Bebe-Ball



Bebe-Ball Manual Suction Device		
Order #	Product Description	B× Qty
	Bebe-Ball with Bebeonkers	100
16-4	Bebe-Ball with Bebeonkers 16-2 (Std)	50
16-5	Bebe-Ball with UP Bebeonkers 16-1 (bub)	50
Be	ebe-Ball with Bebeonkers and Keeper Kove	r
16-4KK	Bebe-Ball with Bebeonkers 16-2 (Std) with Keeper Kover	50
16-5KK	Bebe-Ball with UP Bebeonkers 16-1 (bulb) with Keeper Kover	50



Red ball is our trial color, actual product ion will be blue translucent.

#### Bebe-Ball

- Bulb suction adaptor for Bebeonkers 16-1 and 16-2
- Helps prevent damage that occurs if one-hole device accidently touches the sensitive mucosa of a newborn or preemie
- Disposal after 12 to 24 hours or per hospital policy to prevent infection is recommended
- Use on transports or when power is unavailable and STAT suctioning is required
- Standard bulb syringe has one hole and is inappropriate because it can grab and tear tender tissue too





## Bebe Keeper Kover

Can be used up to 12 hours on the same patient, or per hospital policy.





Bebe Keeper Kover is an effective cover to keep the Bebeonkers clean

Can be used on Bebeonkers I and Bebeonkers II with suction on or off, protecting the device while not in use.




Can also be used with the Bebe-Ball.

Note new blue translucent color

## Parent Guides

- Written for parents in simple and easy to understand language
- NIDCAP inspired book developed by nurses
- · Knowledge of baby's reactions, positive & negative to the NICU
- · Empowers parents to be more involved
- · Assists with discharge teaching for JCAHO requirements
- Parents are signed off as they learn how to care for their baby



Ene	e que promo mengitivo cominioni-
100	AND DESCRIPTION OF THE PARTY OF
	Committee of the control of the cont
	Space Service Control
	Street at the Self-or Secretary
	and the latter property of the party of the latter of the
	NAME OF STREET OFFICE ADDRESS OF THE OWNER, WHEN PERSON NAMED IN
	and the same of th
- 1	total formality
- 3	Annal by Printer many large arrange & P 414
	pat that feelingth collect respective
	nink-ilar
	The state of the s
	The second second second
100	00000000
	the last last last last last last last
0.000	
100	and the second particular in particular
100	
100	alter any or produced part of the part of the
-	to make you be. Margar is divine you represent you got a financial.
	William
	arthr of post
- 4	Mar IC responses
- 46	THE RESIDENCE OF THE PARTY OF T
	Planter III State
-46	Spirit State in Super-representation of State
	MATTER CONTRACTOR
	publication
	Married Color Colo
	Common de la California de Españo
	service to deposit a Property Associate
1.00	
	V 4 1

Order #	Product Description	Bx Qty
	Parent Guide	
17 A	English Guide for Parents with infants in the NICU	50
17AW	English Guide for Parents with infants in the NICU with Hooks	50
178	Spanish Guide for Parents with infants in the NICU	50
178W	Spanish Guide for Parents with infants in the NICU with Hooks	50
17C	French Guide for Parents with infants in the NICU	50
17CW	French Guide for Parents with infants in the NICU with Hooks	5 0



French Available Upon Request

#### .

to in a brokke which will be it 

Registratings:

Begistratings:

Begistrating

As ASEA. There is a place near special wine to which we came yield a similar many to the fine fine in the revenue clock of each page has sign to a lever to be fine unabsensed active care the proof halor. As your first he may be present on the fine section common your near the fine fine fine page of more in the contract of the contra

SECU Staff in health fregunssys.

#### Now it's time to . . .

Get my name ready at house.

Here supplies ready bifupers, duther, butter, mr./
Get a ter next and adjust it and know how to use it.

Pack and a duther for my well halvy care and he that MD know he is duther for my well halvy care and to that MD know he is the staff income if you waste poor bully long commented.



		-	
-			
-			
		=	
		=	

#### Precious Preemie/Newborn Gowns


- · The Precious Preemie/Newborn Gown is Perfect
- One piece gown with 3 Velcro openings to minimize handling
- Opens in the front or back
- Allow easy access to IV lines & electrode wires
- Machine washable cotton/poly blend, flame resistant



" As clinicians on a mission, seeing nurses cutting products to make them work is all the motivation we need to create something better."

-Kenneth Croteau RRT COO

## Wing Sling



- Soft padded brushed nylon & foam for the back of baby's neck
- Flexible stretch netting to create a sling to protect and cradle the affected arm
- Allows for natural movement while holding in a normal secure position







Order #	Pkg Qty			
Wing 5ling				
11A	Arm Immobilizer for broken clavicles	6		

## Privacy Poncho



#### Created for busy units

- · Eliminates screens and their problems
- Provides privacy for Skin2Skin/Kangaroo care
- · Soft fleece material keeps baby warm
- 100% polyester, machine washable, flame resistant
- Large neck opening so nurse can evaluate infant positioning, lines and mom's breast feeding
- Strong Velcro attachments on both sides to hold IV lines and ventilator tubing
- Designed to fit any Mom or Dad

Orden #	Product Description	Cs Qty
	Privacy Poncho	
12A	Privacy cover for breast feeding in the NICU	10
12AO	Privacy cover for breast feeding at home (available online)	10



#### BeBe Chix



Under Pads for the Wee Ones



Order #	Product Description	Pkg Qty
	BeBe Chix Under Pads	
19A	Infant UnderPads 12"×10"	20

These under pads are designed for the little messes from feedings, IV, blood draw drips or just easy clean up. Can be used to protect Heavenly Foam, Cuddle Up, Nest Sets and more. Save on laundering. Inexpensive and very handy.



## \*\*Introducing the Birth N Bath Bag



#### **Gear Birth N Bath Bag**

- a bubbled layer of polyethylene which creates a greenhouse effect to keep amniotic fluid around infant and reduce evaporation and insensible fluid loss.
- adjustable hood/bag fits infant's head/body keeps room air and drafts away to minimize heat loss.
- a pocket on the back to insert a heel warmer at the shoulders for added heat and positions with the head with neck

Check out on next page

## Birth n' Bath Wraps

#### Disposable insulated wrap for wet infants

- Use at birth and/or at bath time
- · Warm & Insulate against insensible water lass and drafts
- To dry, protect & swaddle infant
- · Safe, no plastic inside to stick to wet skin
- · Quilted, soft, fluffy, absorbent & insulated inner liner
- · Plastic waterproof outer layer
- Disposable and economical



- Opened wrap on warmer
- · Place infant in center
- Pull net bandage over head to close and secure

Order #	Product Description	Pkg Qty
	Birth n' Bath Wraps	
0.68nB	Absorbent & Insulating Wraps at birth n bath	20
06BnB UP	Absorbent & Insulating Wraps at birth n bath	20
068nB CL	Clear & Insulating Wrops at birth n bath	20
0.68nB CL	UP Clear & Insulating Wrops at birth n bath	20

## Cuddle Buns


What makes our diapers so special?



The answer .....
The narrow crotch.
Our pads are figure eight in shape.
This helps keep baby positioned to
prevent abduction.

\*\*Ref. article: NYU Child Study Canter Premature Infants: A Later Look

"C orrectly sized diapers is one such simple invention. For example, it was found that preemies who are put in diapers that are too large must keep their legs spread apart in order to accommodate the extra fabric. This resulted in incorrect turn out of the hips and thus later motor problems."





Order #	Product Description	Bx Qty
	Cuddle Buns Developmental Diapers	
21A	Preemie Size fit infants 2-4 lbs or 1300-2000 gms	360
218	Large Preemie fit infants 4-6 lbs or 1850-2750 gms	360
210	Ultra preemie fit infants under 2 lbs or 1300 gms	360

## Cuddle Buns



- A crotch small enough to prevent abduction
- Double elastic that really works to prevent leaks
- A one way liner to keep baby dry





Note: Diaper is too big. Legs are spread wide.

- No gel to swell
- You can easily get a urine specimen
- A cut out to keep diaper away from umbilical lines



"The making of a preemie diaper was very difficult but the dangers to the baby with legs abducted and plastic tearing their sensitive skin kept me on course to find a solution. It was a 5 year struggle to produce a better diaper."


—M. Sharon Rogone RN

## Beary Small

#### Bili-Buns Phototherapy Diapers

The Beary Small Economical Diapers 22A-22B-22C are designed with a narrow crotch and made with SAP for those situations where extra absorbency is required. Made with hook and loop attachments and newly improved all white opaque outer layer and no tape to stick to sensitive skin.

And now our Beary Small Diapers are permeable by photo therapy lights.





Order #	Product Description	8x Qty
	Beary Small Economical Photo Therapy Diapers	
22A	Preemie fit infants 1.5-2.5 bs or 800-1300 gms	360
228	Lg. Preemie fit infants 2.5-4 lbs or 1300-1600 gms	360
22 <i>C</i>	Ultra preemie fit infants under 1.5 lbs or 800 gms	360







## Beary Small

#### Bili-Buns Phototherapy Diapers

Beary Small Bili-Buns

- Still the same but now better
- Bili-Light penetrates diaper
- Now white opaque outer layer
- Hook & loop closure & no tape
- Economical
- Very absorbent with gel
- 3 preemie sizes
- Narrow crotch

Test showing light penetrating the Bili-Buns Diaper through the fluff.





## Neck Rolls

Undershirts no more!!!

Our soft gel-like beads

are flexible & conform

to fit comfortably & allow
neck to remain in a slightly
extended position for

easier ventilation.





Super soft gel beads





rroduct Description	Co OdiA
Bead Neck Rolls Infant Positioning Device	34117
Presmis Neck Roll with Fleece Cover	10
Preemie Neck Roll with Knit Cover	10
Late Term Preemie Neck Roll with Fleece Cover	10
Late Term Preemie Neck Roll with Knit Cover	10
Ultra Preemie Neck Roll with Fleece Cover	10
Ultra Preemie Neck Roll with Knit Cover	10
Newborn Neck Roll with Fleece Cover	10
Newborn Neck Roll with Knit Cover	10
Foam Neck Rolls Infant Positioning Device	
Presmis Neck Roll with Fleece Cover	10
Preemie Neck Roll with Knit Cover	10
Late Term Preemie Neck Roll with Fleece Cover	10
Late Term Preemie Neck Roll with Knit Cover	10
Ultra Preemie Neck Roll with Fleece Cover	10
Ultra Preemie Neck Roll with Knit Cover	10
Newborn Neck Roll with Fleece Cover	10
Newborn Neck Roll with Knit Cover	10
	Preemie Neck Roll with Fleece Cover Preemie Neck Roll with Knit Cover Late Term Preemie Neck Roll with Fleece Cover Late Term Preemie Neck Roll with Fleece Cover Ultra Preemie Neck Roll with Fleece Cover Ultra Preemie Neck Roll with Knit Cover Newborn Neck Roll with Fleece Cover Newborn Neck Roll with Fleece Cover Newborn Neck Roll with Knit Cover Foam Neck Rolls Infant Positioning Device Preemie Neck Roll with Fleece Cover Late Term Preemie Neck Roll with Fleece Cover Ultra Preemie Neck Roll with Fleece Cover Ultra Preemie Neck Roll with Fleece Cover Ultra Preemie Neck Roll with Fleece Cover Newborn Neck Roll with Fleece Cover

## Neck Roll Covers





Neck Rolls with vinyl wipe clean surface & very soft gel-like beads take the place of those rolled-up undershirts, etc. Available in 4 sizes to meet the needs of large to micro-size babies. A fleece cover is included, knit & disposable covers are also available. The knit & fleece covers are machine washable. Place Neck Roll into a cover and place behind infants neck in supine position, be careful not to hyper extend or impair circulation in the neck. These can be used as hand rolls, as well,



Disposabla covar



		Knit Cover
3	38/	3

Knit	Cover

Prints may vary.

Order #	Product Description	Cs Qty
	Neck Rolls Covers Only	
10NE	Preemie Neck Roll Fleece Cover	10
10NEK	Preemie Neck Rolls Knit Covers	10
10NLTE	Late Term Preemie Neck Rolls Fleece Covers	10
10NLTEK	Late Term Preemie Neck Rolls Knit Covers	10
10NNE	Ultra Preemie Neck Rolls Fleece Covers	10
10NNEK	Ultra Preemie Neck Rolls Knit Covers	10
10NBE	Newborn Neck Rolls Fleece Covers	10
10NBEK	Newborn Neck Rolls Knit Covers	10
	Disposable Covers Only	Pkg Qty
10NED	Preemie Neck Rolls Disposable Covers	20
10NLTED	Late Term Preemie Neck Rolls Disposable Covers	20
10NNED	Ultra Preemie Neck Rolls Disposable Covers	20
10NBED	Newborn Neck Rolls Disposoble Covers	20

Notes	
	Ì
	j



## Knowledge is Power

Contains a wealth of knowledge concerning developmental care and outcomes, also safe practices and standards of care for nurses; as well as the cost savings for the NICUs and hospitals.

This document is available to all distributors and also to nurses and NICUs.

Contact Us for More Information (800)676-0462

